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CE

## Preface

Please read the user manual carefully before taking the product into use.

- ✧ This manual contains operations, assembly methods, and simple faults solutions.
- ✧ This manual applies to our model: 0815013
- ✧ This manual contains wheelchair maintenance and self-checking methods, please put it in proper place.
- ✧ Please provide this manual for reference when other people are going to use this wheelchair.
- ✧ The annotations and illustrations in this manual might be slightly different with the real parts due to quality improvement or changing design. Please in kind prevail.
- ✧ Contact with your dealer if there is any ambiguity or question.
- ✧ Improper use of any vehicle may lead to injury. Unsafe driving could harm yourselves and others.
- ✧ The electric wheelchair is intended to comfortably transport persons with walking difficulties or no walking abilities.
- ✧ This electric wheelchair is designed to transport 1 person only.

Symbols used in this manual.

Warning symbol

Follow the instructions next to this symbol closely.

Not paying careful attention to these instructions could result in physical injury or damage to the wheelchair or the environment.

# 1. Product Introduction

## 1.2 Product Images

REF:0815013

### 1.3 Product Information

0815013	SPECIFICATIONS
<b>Frame</b>	Aluminum, Foldable
<b>Controller</b>	Right or Left Adjustable Controller, Dynamic, DC 60A
<b>Motor</b>	450W x 24 V x 2
<b>Battery</b>	2 pcs x 12V 50Ah , lead-acid
<b>Charger</b>	24V 8A
<b>Max Speed</b>	8-10 Km/h
<b>Autonomy*</b>	30-35 km
<b>Gradeability</b>	10°
<b>Front Wheel</b>	9", pneumatic
<b>Rear Wheel</b>	14", pneumatic
<b>Armrest</b>	Detachable & Adjustable in width (48-51cm) & height (16-24cm)
<b>Footrest</b>	Aluminum, Detachable, rotating, adjustable in height & length
<b>Seat</b>	Removable nylon mesh fabric cushion with memory foam lining
<b>Brakes</b>	Electromagnetic Brakes
<b>Other features</b>	Mudguard, seat belt , 3" anti-tilt casters, Adjustable Head support , Adjustable Backrest, Electric seat Tilt , horn, front and rear lights (turning and braking), front and rear lights (flash) with alarm function

\*approximative autonomy as also determined by the weight of occupant, the battery condition, the ambient temperature, the way of driving, the terrain and the wind speed

										
42.5-57.5cm	40-50cm	Adjustable depth in 3 positions 45cm 47.5cm 50cm	60cm	60-69cm	79-105.5 cm	127cm	50cm	14" 9"	150kg	60kg (without batteries), 87kg (with batteries)

### **Product Performance**

This series of electric wheelchairs are powered by battery, driven by motor. Users control direction and adjust speed by joystick controller. It is suitable for application of low speed, good road condition and small slope.

### **Main Structure**

This series of electric wheelchairs are mainly composed of two parts:

- 1) Frame Part: frame, wheel, footrest, armrest, backrest etc.
- 2) Electrical Part: controller, motor, battery.

### **Application Range**

Our electric wheelchair is for disabled and elderly people using as walking vehicle, is suitable for outdoor middle distance use. Do not run on the motorways.

## **2. Safety Instructions**

### **2.1 Driving Notice**

#### **General Driving Cautions**

- ♦ Please keep your hands on the armrest to control the wheelchair.
- ♦ Please practice in parks or other safe open places until you can use the wheelchair skillfully.
- ♦ Fully practice driving in a safe place, to master principle of moving forward, stopping and turning circles.
- ♦ Before driving to the road, please be accompanied by caregivers and confirm it is safe.
- ♦ Please strictly follow traffic rules of the pedestrian, do not regard yourself as vehicle driver.
- ♦ Please drive on the sidewalk and the zebra crossing. Do not drive on vehicle lane.
- ♦ Steadily drive, to avoid driving in 'Z' line or sharp turn.
- ♦ Please keep pressure of pneumatic tire properly. The abnormal tire pressure may cause unsteady driving or excessive lost of current.

#### **Accompanied by Caregivers or Avoid Driving in These Conditions**

- ♦ Drive in bad weather, such as rainy day, heavy fog, strong wind, snow.
- ♦ In case of wheelchair is wet out, wipe the water immediately.
- ♦ Drive on the bad road condition, such as muddy, trail, sand road, gravel, etc.
- ♦ Drive on crowded road.
- ♦ Drive on no fence side ditch, pond, etc.
- ♦ Across the railway.
- ♦ When you have to across railway, pause at "turn off" to confirm it is safe, and make sure that the tires will not be stuck by railway.
- ♦ The electric wheelchair is only for personal use, do not carry people or goods, and do not for a traction purposes.

#### **Precautions for Uphill and Downhill**

- ♦ Avoid driving on steep slopes. The tilt range should be up to 10 degrees.
- ♦ Keep moving forward when going uphill or downhill.
- ♦ Keep speed low when going downhill.
- ♦ Driving on steps is prohibited.
- ♦ Avoid potholes.

#### **Precautions for the Caregivers**

- ♦ Caregivers should confirm that the user's feet are on the footrests proper position and ensure clothes do not attach to the wheels.
- ♦ Caregivers should push the wheelchair to move forward to keep safe on steep slope or long slope.

## 2.2 Other Notice

### Repair and Reform

If the repair and modification is necessary, please contact your dealer. Do not modify by yourself, it may lead to accident or malfunction.

### Protect Environment

In order to protect environment, do not casually discard the abandoned wheelchair and scrapped battery. Please contact your local dealer for replacement of new battery.

### Keep the Storage Environment Dry

Do not storage wheelchair in damp condition, such as near-by bathroom. If wheelchair got wet in the rain, wipe the water immediately.

*Do not replace parts or material randomly.*

*Do not add weight randomly, in case to cause imbalance.*


*Do not use other vehicle to pull or push the wheelchair when somebody sits on the wheelchair or the wheelchair is in manual mode. When there is abnormal sound, contact service immediately.*

### Protection from Electromagnetic Interference

The wheelchair should away from strong magnetic fields and large inductive electrical equipment. Pay attention to that if there are sources of electromagnetic interference nearby, as far as away from those sources to avoid electromagnetic interference

## 3 Adjustment of Wheelchair

- Open the box.
- Remove the back cushion, seat cushion, headrest, accessory box.
- Open the wheelchair
- Adjust the backrest

- Adjust the back straps
- Install the back cushion and seat cushion.

- Install the controller in the desired location.

- Open the accessories box and place the headrest tube on the frame. Adjust the position of the headrest.

- Install the footrests.

**Folding the wheelchair:**

Step one: Remove the battery


Step two: Lift the cushion firmly with your hands and fold the wheelchair

### Adjusting and setting the controller

Adjust the controller to your ideal position by following these steps:

- ✓ Loosen the screw (A) on the armrest support tube.
- ✓ Adjust the controller to the appropriate position.
- ✓ Tighten the screws.
- ✓ Then you can adjust the angle of the controller - move the controller outwards in a circular path. At the same time, you can loosen or tighten the screw to adjust the resistance.

### Assembling and Disassembling the Footrest


Place the footrest in the appropriate place on the wheelchair frame and move it from the outside to the inside. To disassemble it, press the plastic part outwards and lift the footrest.

#### Height adjustment

1. Loosen the screw (A) on the footrest support tube.
2. Adjust the footrest to the appropriate height and tighten the screw.

#### Footrest angle adjustment

1. Turn the footrest upwards.
2. Loosen the lock nut (B).
3. Turn the screw up or down until the pedal reaches the appropriate angle.
4. Tighten the lock nut.

### Adjusting and adjusting the anti-tilt casters


Adjust the angle of the casters by unscrewing the nut.

### Seat belt

Be sure to wear a seat belt. There is an auxiliary guide rail on both sides of the seat frame where you can place the belt.

1. Fasten the belt to the upper groove of the guide rail.
2. Make sure the belt buckle is properly locked.
3. The steps for unlocking the belt are shown in the order ①② in the figure below.

### Switching from manual to electric mode


The product is equipped with levers for switching between manual and electric mode. When the lever is released, the brake will be separated from the motor. In manual mode, the product cannot be controlled electronically.

#### Manual mode:

1. Make sure the power wheelchair is on a flat surface.
2. Turn it off.
3. Turn the levers in.


**Electric mode:**

1. Turn the power wheelchair off.
2. Turn the levers out.

**Backrest tilt adjustment**

The backrest tilt adjustment is controlled by the mechanism indicated by the arrow, by pulling outwards the strap located behind the backrest and it can be adjusted to 90°/100°/110°/120°. When the wheelchair needs to be folded, the backrest can also be folded forward to fit tightly against the seat cushion, thus reducing the volume of the wheelchair.

### Electric Seat Tilt Adjustment


The seat tilt adjustment (0°-25°) will tilt the entire seat of the wheelchair, including the backrest, seat and footrests. When the seat is tilted, the maximum speed of the wheelchair will automatically be reduced or the wheelchair's driving ability will be completely limited.

### Seat width and depth adjustment

You can adjust the seat depth in 3 positions: 45/47.5/50cm using the mechanism indicated by the red arrows below :

You can adjust the seat width 40-50cm using the mechanism below. When adjusting the seat width, first loosen the screws under the seat surface, and then pull the seat plate on both sides, as shown in the picture below.

## 4 Battery and Charger

### 4.1 How to fix the battery box.


### BATTERY INSTALLATION

PLACE THE BATTERY IN TO THE TRAY


TIE THE VELCRO AND PLACE THE BATTERY IN THE BATTERY BOX

ADJUST THE BATTERY TRAY TO THE CORRECT POSITION

CONNECT THE POWER PLUG

REATTACH THE VELCRO AND WRAP THE POWER CORD

THE BATTERY WILL BE AS IN THE PICTURE ABOVE

INSTALL THE BATTERY BOX COVER (BOTTOM)

INSTALL THE BATTERY BOX COVER (TOP)


## **BATTERY REMOVAL**

**REMOVE THE TOP BATTERY COVER BY GRASPING ITS TOP EDGE AND PRESSING DOWNWARDS**

**REMOVE THE BOTTOM BATTERY COVER BY GRASPING ITS TOP EDGE AND PRESSING DOWNWARDS**


**REMOVE THE VELCRO**

**DISCONNECT THE POWER PLUG**

**FIX YOUR FINGERS AROUND THE RING AND WITH YOUR OTHER HAND, GRAB THE VELCRO AND PULL IT OUTWARDS**

**HOLD THE STRAP WITH YOUR HAND**

**REMOVE THE BATTERY**

## 4.2 Charger Requirements


The charger is used for battery charging. During electric wheelchair is charging, do not use it.

### Charger Technical Data

Input Voltage: AC100~240V    Output: 24V/DC    Output Current: 8A

## 4.3 Usage of Charger

You can charge your wheelchair via the control panel. The charging socket is located on the bottom of the control panel.

In order to charge the battery, connect the charger with power supply and battery box's plug. Do as following instructions to complete charging process:

*Step 1: Make sure the charger groove is not blocked.*

*Step 2: Make sure the electric wheelchair is power off.*

*Step 3: Unplug the plugs that connect the battery box and controller.*

*Step 4: Connect the charger's output plug to the battery box's power plug.*

*Step 5: Connect the charger's main plug to the power supply and the red light will be alight. Fully charging needs 8 to 10 hours, do not overcharge more than 24 hours.*

*Step 6: As to prevent cut short battery life, please charge the battery at least once per month when the wheelchair is not being used.*

*Do not stop charging until charge process is finished. Repeat using the battery, which is not fully charged, will shorten the battery life, so the battery should be fully charged as much as possible. When battery is fully charged, the power indicator will turn to green. Do not stop charging before fully charged.*

*After finish charging, turn off the power supply, otherwise the battery will discharge slowly. Do not charge more than 24 hours. Over-charging is dangerous.*

*Users should follow the following rules to avoid charging dangers:*

*The electric wheelchair does not include the charger, please use national standard charger that the output voltage is 24V, 8A.*

*It should be well ventilated when charging. Do not expose the wheelchair to sunlight and humid environment.*

The charging environment temperature range is from 10°C to 50°C. If it is out of the environment temperature range, the battery is unable to function well, and can easily make battery damaged. It is normal for the fan making a sound during charging. It is for cooling the charger, please do not worry about it.

Prevent liquid go into the charger during charging. Moreover, do not place charger on the flammable items, such as fuel, footrest or seat cushion.

## 5 The Usage and Maintenance of Battery

Wrong operation of replacing battery may cause danger of explosion. Only the same or recommended type of battery is suitable for replacement. In addition, please make sure the battery poles are correct. Key points for prolong the battery life: charge frequently, to keep battery power full. It is better to fully charge the battery if the wheelchair is not being used. If stop using for a long time, it is better to charge twice per month.

## 6. Controller Controller Panel

## 7. Maintenance

Maintenance of wheelchair includes clean wheelchair, check wheel and battery, and charge battery. For further maintenance, please contact with your dealer. Our suggestion is checking the wheelchair biannual every year.

## 7.1 Clean and Battery Maintenance

### Clean Wheelchair

Please clean the wheelchair regularly.

Clean the parts, which are frequently touched with user's body (such as seat cushion, armrest, controller) by a clean and slight wet cloth. Do not use organic solvents to clean.

If it is a patient user, the wheelchair should be cleaned once per week. If the wheelchair is used by an infectious victim, it should be cleaned and disinfected by disinfectants.

### Wheels

Check the tires' air pressure and their wear condition regularly. When the tire tread pattern depth reduces to 1mm, please replace new tire.

### Battery

To make sure that the battery is fully charged frequently. In order to prolong battery life, we strongly recommend users do not to charge until battery run out of power total.

## 7.2 Wheelchair Malfunction and Checking

When malfunctions happen during the wheelchair is working, please turn off the power before check.

Symptom: completely lose power, and all the LED lights on the controller panel are off.

### Check Step:

*Step 1:* To check if the controller's plug is loose.

*Step 2:* To check if the connection of controller's plug and battery box is loose. Please reinsert the plug connector (hold the plug when pull out the plug. Do not pull the wire to avoid unnecessary damage to the wire line). After above checking, if the wheelchair is still not able to restore the power, or if users have any question for above checking, please connect with your dealer.

## 7.3 Maintenance Checking

The following is a checking list. Electric wheelchair should be ordinal checked according to our suggestion. When you get on or get off the wheelchair, some self checking is done automatically.

For more your attention, we particularly list these self-checking items in A area.

<b>A area</b>	Before use, please check if following parts are correct: 1.Backrest 2. Armrest 3. Controller position 4.Footrest 5. battery power 6.Clutch/adjusting lever for conversion between manual mode and electric mode
<b>B area</b>	Check the following parts monthly, to avoid original parts losing or wearing 1. Screws 2. Brakes 3. Clutch/adjusting lever for conversion between manual mode and electric mode 4. The front and rear wheels and their tread pattern depth 5. Connectors' of controller and charger
<b>C area</b>	For safety, semiannual overall maintenance is needed.

## 8. Condition for Transport and Storage

Transportation should be avoided moisture and sunlight and away from heat resource.

- ♦ Storage condition:
  - Environment temperature -40°C ~ +55°C;
  - Relative humidity ≤ 80%;
  - Air Pressure 86kPa ~ 106kPa.

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## 10. Warranty

The product is guaranteed for a period of two\* (2) years from the purchasing date.

The warranty covers factory defects & **DOES NOT** cover damages resulting from misuse, poor maintenance, modification, overuse or non-compliance with the product's operating instructions. Also **DOES NOT** cover the parts that can wear out during use or over time like brakes, wheels and plastic parts.

Damages or defects caused by: natural disasters, unauthorized maintenance or repair, supply problems, carriers are not covered by the warranty.

**\*SIX (6) months for the battery**

No service or spare part is covered unless it has been inspected by an authorized service department or the manufacturer of the product. Shipping costs for goods of warranty are to be paid by the user- customer or by dealer

### Warranty Claim Form

<b>BUYERS DETAILS:</b>			
<b>FULL NAME:</b>			
<b>PURCHASING DATE:</b>		<b>PHONE::</b>	
<b>SERIAL NUMBER:</b>			
<b>LOT :</b>			
<b>DEALER DETAILS</b>			
<b>FULL NAME:</b>			
<b>ADDRESS:</b>			
<b>PHONE:</b>			
<b>SIGN AND SIGNATURE</b>			

### INCIDENT REPORTING:

Any serious incident that occurs in relation to use of the device must be reported to the manufacturer and a competent authority of the Member State in which the user and / or patient resides.

### DECLARATION OF CONFORMITY:

We are solely responsible for declaring that the Medical Devices mentioned in this statement are of Low-Risk Class (Class I) and comply with the requirements of the European Regulation 2017/745 and where appropriate, the standards and legislation referred to.